

Please note:

A location = Keizerzaal & Foyer
 B location = Picke 1 room
 C location = Van Berckel Lounge


- = BOOTH AVAILABLE (8m2)
- = BOOTH AVAILABLE (6m2)
- = BOOTH RESERVED
- = BOOTH ON OPTION

Please note:

6m2 booths are numbered:

- 1 – 4
- 10 – 11
- 14 – 21
- 38 – 43

See room measurements below



- Booth 1: Finbiosoft
- Booth 2: Labhelp
- Booth 3: QIAGEN
- Booth 4: Seegene
- Booth 5: Luminex
- Booth 6: GenePOC
- Booth 7: BioMérieux
- Booth 8: Cepheid
- Booth 9: BD Diagnostics
- Booth 10: Mediphos
- Booth 11: R-Biopharm
- Booth 12: Roche royal sponsor
- Booth 13: Roche royal sponsor
- Booth 14: Genetic Signatures
- Booth 15: Curetis
- Booth 16: CYCLERtest & BIOplastics
- Booth 17: Medical Wire
- Booth 18: Perkin Elmer
- Booth 19: Fast track Diagnostics
- Booth 20: Molzym
- Booth 21: Elitech
- Booth 22: PathoFinder B.V.
- Booth 23: Biologio
- Booth 24: Greiner Bio-One
- Booth 25: Multipicom
- Booth 26:
- Booth 27:
- Booth 28:
- Booth 29:
- Booth 30: Qnostics
- Booth 31: Altona Diagnostics
- Booth 32: Abbott B.V.
- Booth 33:
- Booth 34: GenMark Diagnostics (option)
- Booth 35: Mobidiag Ltd
- Booth 36: Siemens Healthcare
- Booth 37: UgenTec
- Booth 38: Hologic
- Booth 39: Copan group
- Booth 40: Speedx
- Booth 41: Diasorin
- Booth 42: Biocartis
- Booth 43: Diagenode Diagnostics

Keizerzaal measurements:

Square Meters:	400
General room measurements:	Width 21,23 / Length 19,90
Measurements lower section:	Width 21,23 / Length 7,60
Ceiling height higher part:	3,85
Ceiling height lower part:	2,70
Window surface:	corner to corner 5,58
Projection screen in the room:	Height 2,97 / Width 5,88
Door Height:	2,22 / Width 2,02

Pillar Keizerzaal to cabinet entrance kitchen: 2,01

Pillar Keizerzaal to wall Bar Maritime: 2,36

Picke 1 measurements:

Square Meters Pické 1:	374
Pické 1 L x W x H:	17 x 23 x 5
Ceiling heights:	6,48
Up to Chandelier:	4,67
Ceiling height rails dividing walls:	5,46
Ceiling height outer ring (wood):	3,38
Ceiling height outer ring:	3,52
Door Pické 1: Height 2,22 / Width 1,95 (between rails 1,72	
Walls landside per surface:	Width 7,42
Walls curves land side:	Width 6,35
Walls middle Pické 1:	Width 6,90
Walls with tapestry:	Width 4,69
Hanging points Max weight:	800 kg per hanging point
Diameter eye hanging point:	0.4
Measurements trusses:	Pické 1: 7,80 x 5,40

Measurements opening translator cabins:

Height 1,40 / Width 1,56 beam in the centre 0,68